

Share With:

**Application:**

P-1 laboratory grinding machine/Metallographic Specimen grinding Polishing Machine is suitable for the polishing of the specimen which has been grinded. The specimen surface is very smooth after processing and can be used to observe and measure the metallographic structure of specimen under microscope. This machine has features of smooth running, low noise, easy operation and maintenance. It can meet polishing requirements of multiple materials. It is widely used in metallographic experiments of mining and industrial enterprises, colleges and universities, scientific research institutions, etc.

**Main specifications:**

Model	P-1
Polishing disc diameter	203mm
Rotating speed	1400r/min (special order)
power supply	AC220v,50Hz
Net weight	15kg
Motor	180w
Dimension	360×450×320mm